# Sistema total de rodilla Triathlon Tritanium<sup>®</sup>

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componentes no cementados.<sup>50</sup>

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- Jackson J.D., Pagnano M.W. (2012) Cement fixation for total knee arthroplasty. En: *The Knee Joint*. Springer, Paris. p. 759-764.
- Dalury, DF. Cementless total knee arthroplasty. *Bone Joint J* 2016;98-B:867-73.
- Webb, JCJ and Spencer, RF. The role of polymethylmethacrylate bone cement in modern orthopaedic surgery. *BJJ. VOL. 89-B, N.º 7, JULIO 2007*
- Barrack, RL. Early Failure of Modern Cemented Stems. *The Journal Of Arthroplasty*. Vol. 15, N.º 8 de ago. de 2000.
- Lettich, T. et al., Primary Total Hip Arthroplasty with an Uncemented Femoral Component: Two- to Seven-Year Results. *The Journal of Arthroplasty* Vol. 22 N.º 7 Supl. 3 2007
- Fleischman, AN et al. Reduced Incidence of Intraoperative Femur Fracture With a Second-Generation Tapered Wedge Stem. *The Journal of Arthroplasty* 32 (2017) 3457-3461
- Registro Nacional de Artroplastias de Reino Unido. XV Informe anual, 2018. <http://www.njrreports.org.uk/Portals/0/PDFdownloads/NJR%2011th%20Annual%20Report%202014.pdf> (Última visita, 11 de nov. de 2018).
- Asociación Australiana de Ortopedia Registro australiano Informe anual, 2018. (Última visita, 11 de nov. de 2018).
- Registro de Artroplastias de los EE. UU.: Informe anual, 2017.
- Memtsoudis SG, Della Valle AG, Besculides MC, Gaber L, Laskin R. (2009) Trends in Demographics, Comorbidity Profiles, InHospital Complications and Mortality Associated with Primary Knee Arthroplasty. *Journal of Arthroplasty*, Jun; 24(4):51827
- Schreurs, BW, Hannik, G. Total Joint Arthroplasty in Younger Patients: Heading for Trouble? *The Lancet* (2017), volumen 389, Número 10077, P. 1374-1375.
- Nam, D., et al. Perioperative and Early Postoperative Comparison of a Modern Cemented and Cementless Total Knee Arthroplasty of the Same Design. *The Journal of Arthroplasty* 32 (2017) 2151-2155
- Nilsson K G, Henricson A, Norgren B, Dalen T. Uncemented HA-coated implant is the optimum fixation for TKA in the young patient. *Clin Orthop Relat Res* 2006; 448: 129-39.
- Meneghini RM, Hanssen A. Cementless Fixation in Total Knee Arthroplasty: Past, Present, and Future. *J Knee Surg*. 2008; 21:307-314.
- Bayley JC, Scott RD, Ewald FC, Holmes GB Jr. Failure of the metal-backed patellar component after total knee replacement. *J Bone Joint Surg Am*. 1988 Jun; 70(5):668-74.
- Bircher, A. et al. Contact Dermatitis. Allergic complications from orthopaedic joint implants: the role of delayed hypersensitivity to benzoyl peroxide in bone cement. *Contact Dermatitis*. Agosto de 2011, 66, 20-26
- Granchi D., Cenni E., Tigani D., Trisolino G., Baldini N., Giunti A. (2008). Sensitivity to implant materials in patients with total knee arthroplasties. *Science Direct*, 29, 1494-1500. doi:10.1016/j.biomaterials.2007.11.038
- Vega, F. Aseptic loosening of a total knee prosthesis caused by delayed hypersensitivity to bone cement. *Annals of Allergy, Asthma & Immunology*. Volumen 117, Número 1, julio de 2016, páginas 89-91
- Keene, R. et al. Occupational Hazards to the Pregnant Orthopaedic Surgeon. *J Bone Joint Surg Am*. 2011; 93:e141(1-5)
- Pacheco, Karin A. Allergy to Surgical Implants. *The Journal of Allergy and Clinical Immunology: In Practice*, Volumen 3, Número 5, 683 – 695
- Nam, D. et al. Cemented Versus Cementless Total Knee Arthroplasty of the Same Modern Design. *J Bone Joint Surg Am*. 2019;101:1185-92
- Sharpe, K., Robinson, K., Cohen, R., Barnett TM., Rastogi, A., Masini, M. Does Implant Fixation Affect Early Return to Function Following Primary Total Knee Arthroplasty. *ISTA 2018*. Londres, RU.
- Sinicrope BJ, Feher AW, Bhimani SJ, Smith LS, Harwin SE, Yakkanti MR, Malkani AL, Increased Aseptic Failures in Cemented versus Cementless Total Knee Arthroplasty in Morbidly Obese Patients. Minimal 5 Year Follow-up, *The Journal of Arthroplasty* (2018), doi: <https://doi.org/10.1016/j.arth.2018.10.016>.
- Miller, AJ, et al., Results of Cemented vs Cementless Primary Total Knee Arthroplasty Using the Same Implant Design. *J Arthroplasty*. 2018 Apr;33(4): 1089-1093. doi: 10.1016/j.arth.2017.11.048. Epub 2 dic. 2017.
- Cohen et al. «Early Clinical Outcomes of a New Cementless Total Knee Arthroplasty Design. *Orthopedics*. 2018.

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26. Buzhardt, P. et al. «Clinical and radiographic results of a highly porous titanium cementless tibial baseplate in TKA». *Bone Joint J* 99-B.SUPP 3(2017): 56. Web. 13 de mar. 2017.»
27. Harwin, SF. et al. Outcomes of Newer Generation Cementless Total Knee Arthroplasty: Beaded Periapatite-Coated vs Highly Porous Titanium-Coated Implant. *J Arthroplasty*. 2017 Jul;32(7):2156-2160. doi: 10.1016/j.arth.2017.01.044. Epub 3 de feb. de 2017.
28. Bhowmik-Stoker et al. «Clinical Performance of a Novel 3D Printed Cementless Titanium Tibial Baseplate, 2-4-year Follow-up» Congreso Mundial de Artroplastia, 4/19-21 de 2018, Roma.
29. King, J.D., et al. The Hidden Cost of Commercial Antibiotic-Loaded Bone Cement: A Systematic Review of Clinical Results and Cost Implications Following Total Knee Arthroplasty. *The Journal of Arthroplasty* 33 9(2018) 3789-3792.
30. Registro de Artroplastia de Rodilla de Suecia. Informe anual 2017. <http://www.myknee.se/en/>. (Última visita, 11 de nov. de 2018).
31. Sultan, A.A., et al. Routine use of commercial antibiotic-loaded bone cement in primary total joint arthroplasty: a critical analysis of the current evidence. *Ann Transl Med* 2019;7(4):73
32. McCalden et al «Comparison of Outcomes and Survivorship Between Patients of Different Age Groups Following TKA». *The Journal of Arthroplasty*, Churchill Livingstone, 24 de julio de 2013, [www.sciencedirect.com/science/article/pii/S0883540313004774](http://www.sciencedirect.com/science/article/pii/S0883540313004774).
33. Meehan, JP, et al. Younger Age Is Associated with a Higher Risk of Early Periprosthetic Joint Infection and Aseptic Mechanical Failure After Total Knee Arthroplasty. *J Bone Joint Surg Am*. 2014; 96:529-35.
34. Mont MA, Gwam C, Newman JM, Chughtai M, Khlopas A, Ramkumar PN, Harwin SF. Outcomes of a newer-generation cementless total knee arthroplasty design in patients less than 50 years of age. *Ann Transl Med* 2017;5(Suppl 3):S24. doi: 10.21037/atm.2017.08.20
35. Newman, J. et al. Cementless Total Knee Arthroplasty in Patients Older Than 75 Years. *The Journal of Knee Surgery*. Vol. 30 N.º 9/2017, P. 930-934.
36. Martin, JR, et al. Morbid Obesity and Total Knee Arthroplasty: A Growing Problem. *Journal of the American Academy of Orthopaedic Surgeons*: Marzo de 2017 - Volumen 25 - Número 3 - p. 188-194. doi: 10.5435/JAAOS-D-15-00684
37. Derman, PB, Fabricant, PD, David, G. The Role of Overweight and Obesity in Relation to the More Rapid Growth of Total Knee Arthroplasty Volume Compared with Total Hip Arthroplasty Volume. *J Bone Joint Surg Am*. 4 jun. 2014;96(11):922-928. Epub 4 de jun. de 2014.
38. Abdel, MA. Et al. Increased aseptic tibial failures in patients with a BMI  $\geq 35$  and well-aligned total knee arthroplasties. *J Arthroplasty*, 30 (2015), pp. 2181-2184
39. Sharpe, K., Robinson, K., Cohen, R., Barnett TM., Cohen, R., Masini, M. Prospective Assessment of the Impact of Obesity on Early Postoperative Outcomes of Cementless Total Knee Arthroplasty AAHKS 2017
40. Sinicrope BJ, Feher AW, Bhimani SJ, Smith LS, Harwin SF, Yakkanti MR, Malkani AL, Increased Aseptic Failures in Cemented versus Cementless Total Knee Arthroplasty in Morbidly Obese Patients. Minimal 5 Year Follow-up, *The Journal of Arthroplasty* (2018), doi: <https://doi.org/10.1016/j.arth.2018.10.016>.
41. Bhowmik-Stoker et al. «Clinical Performance of a Novel 3D Printed Cementless Titanium Tibial Baseplate, 2-4-year Follow-up» Congreso Mundial de Artroplastia, 4/19-21 de 2018, Roma.
42. Firestein GS. Evolving concepts of rheumatoid arthritis. *Nature*. 2003;423(6937):356-361. doi:10.1038/nature01661
43. Helmick CG, et al. Estimates of the Prevalence of Arthritis and Other Rheumatic Conditions in the United States. Part I. *Arthritis Rheumatism*. 2008;58(1):15-25.
44. Patel, N., et al. Outcomes of Cementless Total Knee Arthroplasty in Patients With Rheumatoid Arthritis. *Othopedics*, vol. 41, n.º 2, 2018, p. 103-106
45. Ryd et al. 1995 Ryd L, Albrektsson B E, Carlsson L, Dansgard F, Herberts P, Lindstrand A, Regner L, Toksvig-Larsen S. Roentgen stereophotogrammetric analysis as a predictor of mechanical loosening of knee prostheses. *J Bone Joint Surg Br* 1995; 77(3): 377-83.

46. Nilsson K et al. Evaluation of Micromotion in Cemented vs Uncemented Knee Arthroplasty in Osteoarthritis and Rheumatoid Arthritis. *Journal of Arthroplasty*. Vol. 6. N.º 3. Septiembre de 1991. 265-278.
47. Bhimji et al. The effect of fixation design on micromotion of cementless tibial baseplates. ORS 2012. Póster 1977.
48. Informe de pruebas de Stryker RD-13-107.
49. Protocolo de pruebas de Stryker 92911.
50. Nevelos, J. et al. Design, Migration and Early Clinical Results of the First Mass Produced 3D Printed Cementless Total Knee Implants. Evidencia científica presentada ante la AAOS. Las Vegas, NV, 12-16 de marzo de 2019.
51. Pijls, BG, Plevier J, Nelissen, R. RSA migration of total knee replacements: A systematic review and meta-analysis. *Acta Orthopaedica* 2018; 89 (3): 320-328
52. Dunbar, MJ et al. Stable migration of peri-apatite-coated uncemented tibial components in a multicentre study. *Bone Joint J* 2017;99-B:1596-1602.
53. Robinson JC, Pozen A, Tseng S, Bozic KJ. Variability in costs associated with total hip and knee replacement implants. *J Bone Joint Surg Am*. 2012 Sep 19;94(18): 1693-8.
54. Macario, A. What does one minute of operating room time cost? *Journal of Clinical Anesthesia* (2010) 22, 233–236
55. Volpin A, Khan O, Haddad FS. Theater Cost Is £16/Minute So What Are You Doing Just Standing There? *J Arthroplasty*. Enero de 2016;31(1):22-6. doi: 10.1016/j.arth.2015.08.008. Epub 14 ago. de 2015
56. Chen, A et al. Cementless versus Cemented TKA: Shorter Operative Time and Similar Complications. Eastern Orthopedic Association 2016. New Orleans, LA.
57. Lawrie, C. et al. The Cost of Implanting a Cemented vs. Cementless Total Knee Arthroplasty. Póster n.º 72. AAHKS 2018. Dallas, TX
58. Hansen, EN, et al. Routine use of antibiotic laden bone cement for primary total knee arthroplasty: impact on infecting microbial patterns and resistance profiles. *J Arthroplasty*. Junio de 2014;29(6):1123-7. doi: 10.1016/j.arth.2013.12.004. Epub 10 dic. de 2013.
59. Inabathula, A. et al. Extended Oral Antibiotic Prophylaxis in High-Risk Patients Substantially Reduces Primary Total Hip and Knee Arthroplasty 90-Day Infection Rate. *J Bone Joint Surg Am*. 19 dic. 2018;100(24):2103-2109.
60. Kallala, R., Anderson, P., Morris, S., Haddad, FS. The cost analysis of cemented versus cementless total hip replacement operations on the NHS. *Bone Joint J* 2013;95-B:874-6.

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Copyright © 2019 Stryker